

all-in-ONE™ Enzymatic Enzyme Detergent Surgical Instrument Cleaner with Conditioners

Highly concentrated, fast acting, and cost effective presoak and cleaner.

MJMedical-Inc. website: www.MJMedical-inc.com

Imagine: replacing multiple products, lowering costs, and improving reprocessing outcomes with the all-in-ONE™ Surgical Instrument Cleaner with Conditioners. The all-in-ONE™ rapidly breaks down and removes all forms of proteinaceous bioburden, cleans and conditions the surface, enhances the "passive layer" protection of surgical stainless steel, lubricates moving parts, and then rinses clean. **ONE** gallon of the all-in-ONE™ highly concentrated Surgical Instrument Cleaners with Conditioners will deliver **512 gallons** (1938 liters) of Surgical Instrument Cleaning with Conditioning power. The Dilution rate is .25 to .5 ounces per gallon (1.92 ml per liter to 3.85 ml per liter) depending on level of bioburden, encrustation and **STAINS**. The **Neutral pH** all-in-ONE™ is non-irritating to eyes, skin, and is **100% biodegradable**.

The all-in-ONE™ easy FOAM-it™ lays down a thick blanket of long lasting foam to hydrate and **remove** all forms of proteinaceous debris, hard water deposits, and encrusted **STAINS**, for fast acting cleaning action. When water is added the 'foam' dissipates providing a clear '**free-rinsing**' cleaning solution.

all-in-ONE™ Surgical Instrument Cleaner with Conditioners

ONE product for Cleaning and Conditioning
surgical instruments and scopes.

ONE application for thoroughly
Cleaning and Conditioning.

Apply it **ONCE**...

It's all you really need for Cleaning and Conditioning
surgical instruments, utensils, and scopes.

Replace Cleaning and Conditioning products for:
removing all forms of proteinaceous bioburden,
pre-cleaning and soaking to remove **STAINS**,
manual cleaning, Ultrasonic cleaning, and
automated Surgical Instrument Washing.



easy FOAM-it™
Surgical Instrument
Cleaner and Conditioners



all-in-ONE™
Surgical Instrument
Cleaner with Conditioners

The all-in-ONE™ Surgical Instruments Cleaner with Conditioners delivers a unique formulation of detergents and a complete (**four** multi-tiered "high-level" **enzymes**: Protease, Amylase, Lipase, Carbohydrase) within a **Neutral pH** surfactant chemical complex. The all-in-ONE™ cleaners will: break down bioburden and clean the surface, remove **STAINS**, enhance the "passive layer" protection of surgical stainless steel, and lubricate moving parts. The **highly concentrated** all-in-ONE™ is aggressive against your most difficult cleaning challenge and guaranteed to be safe for all types of metal, rigid and flexible scopes, lenses, microsurgical instrumentation and reprocessing equipment. The all-in-ONE™ delivers more cleaning with conditioning power than any other Enzymatic Surgical Instrument Cleaner. Compare the cleaning performance of the all-in-ONE™ to the lesser effective 'enzyme' detergents against your most difficult reprocessing challenge. You will discover the superior cleaning performance of the all-in-ONE™ Surgical Instruments Cleaner with Conditioners.

UPS-FS-32 all-in-ONE™ easy FOAM-it™ 12/case
Ergo UPS1-1 all-in-ONE™: 4 gallons/case
Ergo UPS1-5 all-in-ONE™: 1 ea. 5 gallon
Ergo UPS1-15 all-in-ONE™: 1 ea. 15 gallon
Ergo UPS1-30 all-in-ONE™: 1 ea. 30 gallon
Ergo UPS1-55 all-in-ONE™: 1 ea. 55 gallon

Dilution Rate: .25 ounce to .5 ounce per gallon depending on the level of bio-burden and **STAINS**.

MJMedical-Inc.

P.O. Box 78678, Charlotte, NC 28271-7038 email: mjhardie@windstream.net

Office: 704-608-9174 Fax: 704-814-4029 website: www.MJMedical-inc.com

Technical Support: Support@ergo-Logistics.com 1-509-747-5027 Website: ergo-Logistics.com